

Oil Bath



- Exterior construction is made of mild quality steel with powder coating finish in standard model.
- In GMP model Exterior construction comprises of matt finish full stainless steel 304 GRADE high quality is used.
- High density glass wool insulation is used between inner & outer wall during manufacturing process to prevent heat dissipation
- Capillary thermostat is used to control temperature with accuracy of 2 °C.
- High Quality Inversion type of heater is used for heating.
- Instrument is supplied with Lid.
- Temperature range 5° C above ambient to 300° C.
- Operates on 230 Volts AC Single phase 50 Hz.

OPTIONAL ACCESSORIES

Dual display microprocessor based PID temperature controller PT 100 Sensor.

MODEL	FLASK CAPACITY	WORKING SIZE (DIAMETER x HEIGHT) in Inch	WATTS
ARPHI AOB 1K	1 LITER	7" X 7"	1000
ARPHI AOB 2K	2 LITERS	8" X 8"	1000
ARPHI AOB 3K	3 LITERS	9" X 9"	1500
ARPHI AOB 5K	5 LITERS	10" X 10"	2000
ARPHI AOB 10K	10 LITERS	14" X 14"	2000
ARPHI AOB 20K	20 LITERS	16" X 16"	2500

Arphi Electronics Pvt. Ltd
32, Corporate Avaneue, A wing,
4th Floor, Near Paper Box Industries
Off Mahakali Caves Road, Andheri (East)
Mumbai-400 093



info@arphi.com



www.arphi.com



022-49702305

022-49702306

022-49702307