

Water Bath



- Double wall leakproof ARGON welded construction with inner chamber made of 304 high grade quality stainless steel.
 - Exterior construction is made of mild quality steel with powder coating finish in standard model. In GMP model matt finish of fully stainless steel is used.
 - High density glass wool insulation is used between inner & outer wall during manufacturing process to prevent heat dissipation.
 - U shaped immersion type heater is used for heating.
 - Capillary thermostat is used to control temperature with accuracy of 2°C .
 - Temperature range 5°C above ambient to 99°C .
 - Instrument is supplied without Lid.
- Operates on 230 Volts AC Single phase 50 Hz.
- OPTIONAL ACCESSORIES**
- A) Dual display microprocessor based PID temperature controller with temperature accuracy of 1°C or better.
- B) Low level water CUT OFF device with suitable sensor for protection of heater lid cover.

Round

MODEL	FLASK CAPACITY	WORKING SIZE (DIAMETER X HEIGHT)	RATING IN WATTS
ARPHI AWBRo-1K	1 LITER	7" X 7"	1000
ARPHI AWBRo-2K	2 LITERS	8" X 8"	1000
ARPHI AWBRo-3K	3 LITERS	9" X 9"	1500
ARPHI AWBRo-5K	5 LITERS	10" X 10"	2000
ARPHI AWBRo-10K	10 LITERS	14" X 14"	2000
ARPHI AWBRo-20K	20 LITERS	16" X 16"	2500

Rectangular

MODEL	FLASK CAPACITY	WORKING SIZE (DIAMETER X HEIGHT) in Inch	RATING IN WATTS
ARPHI AWBRo-1K	1 LITER	7" X 7"	1000
ARPHI AWBRo-2K	2 LITERS	8" X 8"	1000
ARPHI AWBRo-3K	3 LITERS	9" X 9"	1500
ARPHI AWBRo-5K	5 LITERS	10" X 10"	2000
ARPHI AWBRo-10K	10 LITERS	14" X 14"	2000
ARPHI AWBRo-20K	20 LITERS	16" X 16"	2500

Arphi Electronics Pvt. Ltd
32, Corporate Avenue, A wing,
4th Floor, Near Paper Box Industries
Off Mahakali Caves Road, Andheri (East)
Mumbai-400 093



info@arphi.com



www.arphi.com



022-49702305
022-49702306
022-49702307