

Incubator



- High quality 304 GRADE stainless steel is used for inner construction
- Exterior construction is made of mild quality steel with powder coating finish of fully stainless steel is used.
- High quality Acrylic transparent is used in instrument to have clear view thermal conditions inside the chamber.
- Fully insulated door comprises of Sturdy High-quality stainless-steel hinges, latch handle lock. A full gasket towards door acts as a perfect sealant.
- High quality & grade NICHROME wire are used for heating elements.
- Air ventilators are placed appropriately on the incubators for removal of heat gases & fumes.
- Capillary thermostat is used to control temperature with accuracy of $\pm 2^{\circ}\text{C}$.
- Galvanized wire mesh is used for shelves.
- Temperature range 5°C above ambient to 60°C .
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	SIZE in Inch	CAPACITY (In ltr)	NUMBER of SHELVES
ARPHI LIBh-28	12" X 12" X 12"	28Liters	2
ARPHI LIBh-45	14" X 14" X 14"	45Liters	2
ARPHI LIBh-95	18" X 18" X 18"	95Liters	2
ARPHI LIBh-125	18" X 18" X 24"	125Liters	3
ARPHI LIBh-225	24" X 24" X 24"	225Liters	3
ARPHI LIBh-250	24" X 18" X 36"	250Liters	3
ARPHI LIBh-340	24" X 24" X 36"	340Liters	4
ARPHI LIBh-470	24" X 24" X 50"	470Liters	4

Arphi Electronics Pvt. Ltd
32, Corporate Avenue, A wing,
4th Floor, Near Paper Box Industries
Off Mahakali Caves Road, Andheri (East)
Mumbai-400 093



info@arphi.com



www.arphi.com



022-49702305
022-49702306
022-49702307