

## Vertical Autoclave



- High Quality Stainless Steel used for making inside & Outside chambers.
- Advance microprocessor controller for setting time and temperature .
- upto 122 °C Temperature range .
- Precise monitoring using PT 100 Sensor.
- Alarm Indication on low water level & high temperature.
- Buzzer indication for Cycle end .
- Safety valve automatic exahsut on high pressure .

MODEL	WORKING CHAMBER SIZE in Inch	Heater in Watt	Ltr ( in Capacity)	AUTOCLAVE STANDARD (MANUALLY OPERATED)	AUTOCLAVE SEMI-AUTOMATIC	AUTOCLAVE FULLY AUTOMATIC
ARPHI AUTOCLAVE A05	18 x 18"	2000 KW	22	Available	Available	Available
ARPHI AUTOCLAVE A06	12 x 20 "	2000 KW	35	Available	Available	Available
ARPHI AUTOCLAVE A07	14 x 20 "	3000 KW	52	Available	Available	Available
ARPHI AUTOCLAVE A08	16 x 24 "	4000 KW	90	Available	Available	Available
ARPHI AUTOCLAVE A09	18 x 24 "	4000 KW	95	Available	Available	Available
ARPHI AUTOCLAVE A10	18 x 30 "	4500 KW	120	Available	Available	Available
ARPHI AUTOCLAVE A11	22 x 30 "	5000 KW	175	Available	Available	Available
ARPHI AUTOCLAVE A12	22 x 36 "	6000 KW	200	Available	Available	Available

**Arphi** Electronics Pvt. Ltd  
 32, Corporate Avauue, A wing,  
 4th Floor,Near Paper Box Industries  
 Off Mahakali Caves Road, Andheri (East)  
 Mumbai-400 093



info@arphi.com



www.arphi.com



022-49702305  
 022-49702306  
 022-49702307